

sensicam uv

uv sensitive digital 12bit CCD camera system

- excellent quantum efficiency in UV (QE~10% @ 193nm)
- electron multiplication gain from 1 to 500
- high resolution of 1004 x 1002 pixel
- extremely low noise $< 1e^-$ @ gain > 50
- 12bit dynamic range @ gain = 1
- shutter / exposure times from 75 μ s to 1h
- thermoelectrical cooling down to -12°C



sensicam uv

This uv sensitive cooled digital 12bit camera system comprises advanced CCD and electronics technology. With this on-chip multiplication of the light signal the readout noise of the camera can be neglected ($< 1 e^- \text{ rms @ gain } > 50$). With its excellent resolution of 1004x1002pixel this high performance cooled digital 12 bit CCD camera system is best suited for extreme low light applications. The system features thermo-electrical cooling of the image sensor (down to $-16 \text{ }^\circ\text{C}$) and an outstanding quantum efficiency in the UV range of the spectrum ($\text{QE} \approx 15\% \text{ @ } 254\text{nm}$ and $\text{QE} \approx 10\% \text{ @ } 193\text{nm}$). Exposure time modes (software selectable) range from $75 \mu\text{s} - 1\text{h}$. A high speed serial data link connects the system to the PC (fiber optic link available).

technical data

| | unit | setpoint | sensicam uv |
|-------------------------------------|---------------------------------|---|---------------------------------------|
| resolution (hor x ver) ¹ | pixel | | 1004x1002 |
| pixel size (hor x ver) | μm^2 | | 8.0 x 8.0 |
| sensor format/ diagonal | mm^2 / mm | | 8.03 x 8.02/ 11.35 |
| quantum efficiency | % | @ 610nm typ. @ 254nm typ. @ 193nm appr. | 65 15 10 |
| full well capacity | e^- | | 40 000 |
| image sensor | | | TC285SPD |
| dynamic range | dB | @ CCD + camera | 65 |
| dynamic range A/D ² | bit | | 12 |
| readout noise | $e^- \text{ rms}$ | @em gain = 1 @em gain > 50 | 22 <1 |
| imaging frequency, frame rate | fps | @full frame | 13 |
| pixel scan rate | MHz | | 16 |
| A/D conversion factor | e^-/count | @em gain = 1 | 9.8 |
| spectral range | nm | | 190..1100 |
| exposure time | s | | $75\mu\text{s}..1\text{h}$ |
| anti-blooming factor | | typical | 1000 |
| smear | % | | 0.6 |
| binning horizontal | pixel | | 1, 2, 4, 8 |
| binning vertical | pixel | full resolution for 992pixel ver | 1, 2 1, 2, 4, 8, 16, 32 |
| dark current | $e^-/\text{pixel}\cdot\text{s}$ | @-15 °C typical | 0.9 |
| region of interest | pixel | for 992pixel ver | down to 32 |
| charge multiplication | | 9 steps | 1, 2, 5, 10, 20, 50, 100, 200, 500 |

technical data

| | | | |
|--|--------------------|--------------------------------|--|
| non linearity | % | full temperature range @gain=1 | <1 |
| uniformity darkness DSNU ³ | e ⁻ rms | @ 90% center zone & gain=1 | 2 |
| uniformity brightness PRNU ⁴ | % rms | typical @ 610nm | 0.6 |
| trigger, auxiliary signals | | internal/external | software / TTL level |
| power consumption | W | typical | 36 |
| power supply | VAC | | 90..260 |
| mechanical dimensions camera (w x h x l) | mm ³ | | 93 x 78 x 210 |
| mechanical dimensions power supply (w x h x l) | mm ³ | | 84 x 50 x 155 |
| weight | kg | | 1.6 |
| operating temperature range | °C | | +5..+40 |
| operating humidity range | % | | 10..90 |
| storage temperature range | °C | | -20..+70 |
| optical input | | | c-mount, Nikon f-mount |
| optical input window | | | fused silica |
| data interface | | | PCI local bus, Rev 2.1, burst rate 132 MByte/s |
| CE certified | | | yes |
| cooled CCD ⁵ | °C | | -12 |
| cooling method | | | 2 stage Peltier cooler with forced air cooling |

software:

Camware software for camera control, image acquisition and archiving of images in various file formats, WindowsNT, 2K and XP, 32bit-dynamic link library (DLL) is available for user customisation and integration on PC platforms (software development kit - SDK), software is operational in either single mode or with built-in recorder functions, drivers for popular third party software packages are available (see website)

options:


custom-made versions

- [1] horizontal versus vertical
- [2] Analog-to-Digital-converter
- [3] dark signal non-uniformity
- [4] photo reponse non-uniformity
- [5] default setting, programmable to -20°C

frame rate table [frames per second]

| vertical resolution [pixel] | vertical binning | frame rate [fps] |
|-----------------------------|------------------|------------------|
| 1002 | 1 | 12.9 |
| 501 | 2 | 25.2 |
| 248 | 4 | 47.9 |
| 124 | 8 | 86.6 |
| 62 | 16 | 86.6 |
| 31 | 32 | 108.4 |

quantum efficiency



(TC285SPD qe curves as measured by PCO AG)

areas of application

laser induced fluorescence ■ metrology ■ fluorescence spectroscopy (down to UV) ■ bioluminescence ■ chemoluminescence ■ low light level imaging

contact

PCO AG
Donaupark 11
93309 Kelheim, Germany

fon +49 (0)9441 2005 50
fax +49 (0)9441 2005 20
info@pco.de
www.pco.de

The Cooke Corporation
6930 Metroplex Drive
Romulus, Michigan 48174
USA

fon +1 248 276 8820
fax +1 248 276 8825
info@cookecorp.com
www.cookecorp.com

sensicam uv preliminary 11/2006
subject to changes without prior notice

© PCO AG, Kelheim